

# Dyslexia-associated protein KIAA0319-like protein

Catalog No: #AP84668



Package Size: #AP84668-1 50ug #AP84668-2 100ug #AP84668-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dyslexia-associated protein KIAA0319-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Adeno-associated virus receptor,AAVR,Aavr
Accession No.	Uniprot ID:Q8K135
Uniprot	Q8K135
GeneID	100317;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Kiaa0319!

Note: This product is for in vitro research use only