

Heterogeneous nuclear ribonucleoprotein K

Catalog No: #AP84670



Package Size: #AP84670-1 50ug #AP84670-2 100ug #AP84670-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heterogeneous nuclear ribonucleoprotein K
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	hnRNP K,dC stretch-binding protein,CSBP,Hnrpk
Accession No.	Uniprot ID:P61980
Uniprot	P61980
GeneID	117282;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hnrpk

Note: This product is for in vitro research use only