

## Lysine-specific histone demethylase 1A

Catalog No: #AP84674



Package Size: #AP84674-1 50ug #AP84674-2 100ug #AP84674-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lysine-specific histone demethylase 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	BRAF35-HDAC complex protein BHC110,Flavin-containing amine oxidase domain-containing protein 2,AOF2,KDM1,KIAA0601,LSD1
Accession No.	Uniprot ID:O60341
Uniprot	O60341
GeneID	23028;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:KDM1A

Note: This product is for in vitro research use only