

Keratin, type II cytoskeletal 8

Catalog No: #AP84675



Package Size: #AP84675-1 50ug #AP84675-2 100ug #AP84675-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type II cytoskeletal 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytokeratin endo A,Cytokeratin-8,CK-8,Keratin-8,K8,Type-II keratin Kb8,Krt2-8
Accession No.	Uniprot ID:Q10758
Uniprot	Q10758
GeneID	25626;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Krt8

Note: This product is for in vitro research use only