

6-phosphofructokinase, muscle type

Catalog No: #AP84676



Package Size: #AP84676-1 50ug #AP84676-2 100ug #AP84676-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	6-phosphofructokinase, muscle type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Phosphofructo-1-kinase isozyme A,PFK-A,Phosphofructokinase-M,Phosphofructokinase 1,Phosphohexokinase,Pfk-m
Accession No.	Uniprot ID:P47858
Uniprot	P47858
GeneID	65152;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pfkm

Note: This product is for in vitro research use only