

Immunoglobulin superfamily member 1

Catalog No: #AP84679



Package Size: #AP84679-1 50ug #AP84679-2 100ug #AP84679-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Immunoglobulin superfamily member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IgSF1, Immunoglobulin-like domain-containing protein 1, Inhibin-binding protein, InhBP, Pituitary gland-specific factor 2, p120, IGDC1, KIAA0364, PGSF2
Accession No.	Uniprot ID: Q8N6C5
Uniprot	Q8N6C5
GeneID	3547;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IGSF1

Note: This product is for in vitro research use only