

Interferon-induced transmembrane protein 3

Catalog No: #AP84682



Package Size: #AP84682-1 50ug #AP84682-2 100ug #AP84682-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon-induced transmembrane protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Dispanin subfamily A member 2b,DSPA2b,Interferon-inducible protein 1-8U
Accession No.	Uniprot ID:Q01628
Uniprot	Q01628
GeneID	10410;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IFITM3

Note: This product is for in vitro research use only