

Interferon alpha-14

Catalog No: #AP84683



Package Size: #AP84683-1 50ug #AP84683-2 100ug #AP84683-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon alpha-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IFN-alpha-14, Interferon alpha-H, LeIF H, Interferon lambda-2-H
Accession No.	Uniprot ID:P01570
Uniprot	P01570
GeneID	3448;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IFNA14

Note: This product is for in vitro research use only