

Interferon alpha-2

Catalog No: #AP84684



Package Size: #AP84684-1 50ug #AP84684-2 100ug #AP84684-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon alpha-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IFN-alpha-2, Interferon alpha-A, LeIF A, IFNA2A, IFNA2B, IFNA2C
Accession No.	Uniprot ID:P01563
Uniprot	P01563
GeneID	3440;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IFNA2

Note: This product is for in vitro research use only