

Kynurenine--oxoglutarate transaminase 1

Catalog No: #AP84686



Package Size: #AP84686-1 50ug #AP84686-2 100ug #AP84686-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Kynurenine--oxoglutarate transaminase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cysteine-S-conjugate beta-lyase, Glutamine transaminase K, GTK, Glutamine--phenylpyruvate transaminase, Kynurenine aminotransferase 1, Kynurenine aminotransferase I, KATI, Kynurenine--oxoglutarate transaminase I, Ccbl1, Kat
Accession No.	Uniprot ID:Q8BTY1
Uniprot	Q8BTY1
GeneID	70266;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Kyat1

Note: This product is for in vitro research use only