

Serine protease inhibitor Kazal-type 7

Catalog No: #AP84689



Package Size: #AP84689-1 50ug #AP84689-2 100ug #AP84689-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serine protease inhibitor Kazal-type 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Esophagus cancer-related gene 2 protein, ECRG-2, Ecg2
Accession No.	Uniprot ID:Q6IE51
Uniprot	Q6IE51
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Spink7

Note: This product is for in vitro research use only