

Non-histone chromosomal protein HMG-14

Catalog No: #AP84693



Package Size: #AP84693-1 50ug #AP84693-2 100ug #AP84693-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Non-histone chromosomal protein HMG-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	High mobility group nucleosome-binding domain-containing protein 1,Hmg-14,Hmg14
Accession No.	Uniprot ID:P18608
Uniprot	P18608
GeneID	15312;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hmgn1

Note: This product is for in vitro research use only