

# Tyrosine-protein kinase SYK

Catalog No: #AP84698



Package Size: #AP84698-1 50ug #AP84698-2 100ug #AP84698-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tyrosine-protein kinase SYK
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Spleen tyrosine kinase,p72-Syk
Accession No.	Uniprot ID:P43405
Uniprot	P43405
GeneID	6850;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SYK

Note: This product is for in vitro research use only