

Integrin beta-3

Catalog No: #AP84701



Package Size: #AP84701-1 50ug #AP84701-2 100ug #AP84701-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Integrin beta-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Platelet membrane glycoprotein IIIa, GPIIIa, CD61
Accession No.	Uniprot ID:O54890
Uniprot	O54890
GeneID	16416;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Itgb3

Note: This product is for in vitro research use only