

## Inositol monophosphatase 1

Catalog No: #AP84706



Package Size: #AP84706-1 50ug #AP84706-2 100ug #AP84706-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Inositol monophosphatase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IMP 1,D-galactose 1-phosphate phosphatase,Inositol-1(or 4)-monophosphatase 1,Lithium-sensitive myo-inositol monophosphatase A1,IMPA
Accession No.	Uniprot ID:P29218
Uniprot	P29218
GeneID	3612;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IMPA1

Note: This product is for in vitro research use only