

Thyroxine 5-deiodinase

Catalog No: #AP84708



Package Size: #AP84708-1 50ug #AP84708-2 100ug #AP84708-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thyroxine 5-deiodinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	1.21.99.3,5DIII,DIOIII,Type 3 DI,Type III iodothyronine deiodinase,Itid3,Txd3
Accession No.	Uniprot ID:P49897
Uniprot	P49897
GeneID	29475;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dio3

Note: This product is for in vitro research use only