

Interleukin-36 gamma

Catalog No: #AP84710



Package Size: #AP84710-1 50ug #AP84710-2 100ug #AP84710-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-36 gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Interleukin-1 family member 9,IL-1F9,Il1f9
Accession No.	Uniprot ID:Q8R460
Uniprot	Q8R460
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Il36g

Note: This product is for in vitro research use only