

Leukemia inhibitory factor receptor

Catalog No: #AP84717



Package Size: #AP84717-1 50ug #AP84717-2 100ug #AP84717-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Leukemia inhibitory factor receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	LIF receptor,LIF-R,CD118
Accession No.	Uniprot ID:Q5XNR9
Uniprot	Q5XNR9
GeneID	449478;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LIFR

Note: This product is for in vitro research use only