

Leukocyte immunoglobulin-like receptor subfamily A member 1

Catalog No: #AP84718

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84718-1 50ug #AP84718-2 100ug #AP84718-3 1mg

Description

Product Name	Leukocyte immunoglobulin-like receptor subfamily A member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD85 antigen-like family member I,Leukocyte immunoglobulin-like receptor 6,LIR-6,CD85i,LIR6
Accession No.	Uniprot ID:O75019
Uniprot	O75019
GeneID	11024;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LILRA1

Note: This product is for in vitro research use only