

Lipopolysaccharide-induced tumor necrosis factor-alpha factor homolog

Catalog No: #AP84719

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84719-1 50ug #AP84719-2 100ug #AP84719-3 1mg

Description

Product Name	Lipopolysaccharide-induced tumor necrosis factor-alpha factor homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	LPS-induced TNF-alpha factor homolog,Estrogen-enhanced transcript protein 1,Eet-1
Accession No.	Uniprot ID:P0C0T0
Uniprot	P0C0T0
GeneID	65161;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Litaf

Note: This product is for in vitro research use only