

Transcription factor AP-1

Catalog No: #AP84722



Package Size: #AP84722-1 50ug #AP84722-2 100ug #AP84722-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor AP-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Activator protein 1,AP1,Proto-oncogene c-Jun,V-jun avian sarcoma virus 17 oncogene homolog,p39
Accession No.	Uniprot ID:P05412
Uniprot	P05412
GeneID	3725;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:JUN

Note: This product is for in vitro research use only