

Protein KIBRA

Catalog No: #AP84725



Package Size: #AP84725-1 50ug #AP84725-2 100ug #AP84725-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein KIBRA
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HBeAg-binding protein 3,Kidney and brain protein,KIBRA,WW domain-containing protein 1,KIAA0869
Accession No.	Uniprot ID:Q8IX03
Uniprot	Q8IX03
GeneID	23286;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:WWC1

Note: This product is for in vitro research use only