

Leptin

Catalog No: #AP84730



Package Size: #AP84730-1 50ug #AP84730-2 100ug #AP84730-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Leptin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Obesity factor,OB,OBS
Accession No.	Uniprot ID:Q29406
Uniprot	Q29406
GeneID	396832;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LEP

Note: This product is for in vitro research use only