

Inter-alpha-trypsin inhibitor heavy chain H3

Catalog No: #AP84733



Package Size: #AP84733-1 50ug #AP84733-2 100ug #AP84733-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--------------------------------------------------------|
| Product Name | Inter-alpha-trypsin inhibitor heavy chain H3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | ITI heavy chain H3 |
| Accession No. | Uniprot ID:Q63416 |
| Uniprot | Q63416 |
| GeneID | 50693; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Itih3

Note: This product is for in vitro research use only