

Inositol monophosphatase 1

Catalog No: #AP84735

Package Size: #AP84735-1 50ug #AP84735-2 100ug #AP84735-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inositol monophosphatase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IMP 1,D-galactose 1-phosphate phosphatase,Inositol-1(or 4)-monophosphatase 1,Lithium-sensitive myo-inositol monophosphatase A1,IMPA
Accession No.	Uniprot ID:P20456
Uniprot	P20456
GeneID	281865;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IMPA1

Note: This product is for in vitro research use only