

Mitogen-activated protein kinase kinase kinase kinase

1

Catalog No: #AP84740

Package Size: #AP84740-1 50ug #AP84740-2 100ug #AP84740-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitogen-activated protein kinase kinase kinase kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Hematopoietic progenitor kinase,HPK,MAPK/ERK kinase kinase kinase 1,MEK kinase kinase 1,Hpk1
Accession No.	Uniprot ID:P70218
Uniprot	P70218
GeneID	26411;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Map4k1

Note: This product is for in vitro research use only