

Lysophosphatidic acid receptor 1

Catalog No: #AP84742



Package Size: #AP84742-1 50ug #AP84742-2 100ug #AP84742-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysophosphatidic acid receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LPA receptor 1,Lysophosphatidic acid receptor Edg-2,Rec1.3,VZG-1,Edg2,Gpcr26,Lpa1,Vzg1
Accession No.	Uniprot ID:P61793
Uniprot	P61793
GeneID	14745;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lpar1

Note: This product is for in vitro research use only