

# Low-density lipoprotein receptor-related protein 1B

Catalog No: #AP84743



Package Size: #AP84743-1 50ug #AP84743-2 100ug #AP84743-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                                         |
|--------------------|-----------------------------------------------------------------------------------------|
| Product Name       | Low-density lipoprotein receptor-related protein 1B                                     |
| Brief Description  | Recombinant Protein                                                                     |
| Host Species       | E.coli                                                                                  |
| Purification       | Greater than 90% by SDS-PAGE                                                            |
| Species Reactivity | Ms                                                                                      |
| Other Names        | LRP-1B,Low-density lipoprotein receptor-related protein-deleted in tumor,LRP-DIT,Lrpdit |
| Accession No.      | Uniprot ID:Q9JI18                                                                       |
| Uniprot            | Q9JI18                                                                                  |
| Tag Info           | His                                                                                     |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                                  |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Lrp1b

Note: This product is for in vitro research use only