

Lipocalin-1

Catalog No: #AP84744



Package Size: #AP84744-1 50ug #AP84744-2 100ug #AP84744-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lipocalin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Tear lipocalin, Tlc, Tear prealbumin, TP, Von Ebner gland protein, VEG protein, VEGP
Accession No.	Uniprot ID:P31025
Uniprot	P31025
GeneID	3933;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LCN1

Note: This product is for in vitro research use only