

# Early placenta insulin-like peptide

Catalog No: #AP84750



Package Size: #AP84750-1 50ug #AP84750-2 100ug #AP84750-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Early placenta insulin-like peptide
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EPIL,Insulin-like peptide 4,Placentin
Accession No.	Uniprot ID:Q14641
Uniprot	Q14641
GeneID	3641;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:INSL4

Note: This product is for in vitro research use only