

Integrin alpha-V

Catalog No: #AP84752



Package Size: #AP84752-1 50ug #AP84752-2 100ug #AP84752-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrin alpha-V
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Vitronectin receptor subunit alpha,CD51,MSK8,VNRA
Accession No.	Uniprot ID:P06756
Uniprot	P06756
GeneID	3685;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ITGAV

Note: This product is for in vitro research use only