

## Kallikrein-8

Catalog No: #AP84753



Package Size: #AP84753-1 50ug #AP84753-2 100ug #AP84753-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Kallikrein-8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	mK8,Neuropsin,NP,Serine protease 19,Nrpn,Prss19
Accession No.	Uniprot ID:Q61955
Uniprot	Q61955
GeneID	259277;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Klk8

Note: This product is for in vitro research use only