

Protein kinase C epsilon type

Catalog No: #AP84754



Package Size: #AP84754-1 50ug #AP84754-2 100ug #AP84754-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein kinase C epsilon type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	nPKC-epsilon, Pkce
Accession No.	Uniprot ID:P09216
Uniprot	P09216
GeneID	29340;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Prkce

Note: This product is for in vitro research use only