

Creatine kinase B-type

Catalog No: #AP84755



Package Size: #AP84755-1 50ug #AP84755-2 100ug #AP84755-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Creatine kinase B-type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	B-CK,Creatine kinase B chain
Accession No.	Uniprot ID:P05124
Uniprot	P05124
GeneID	100855552;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CKB

Note: This product is for in vitro research use only