

Inter-alpha-trypsin inhibitor heavy chain H3

Catalog No: #AP84758



Package Size: #AP84758-1 50ug #AP84758-2 100ug #AP84758-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Inter-alpha-trypsin inhibitor heavy chain H3 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | ITI heavy chain H3, Serum-derived hyaluronan-associated protein, SHAP |
| Accession No. | Uniprot ID:Q06033 |
| Uniprot | Q06033 |
| GeneID | 3699; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ITIH3

Note: This product is for in vitro research use only