

# Serine/threonine-protein kinase LMTK3

Catalog No: #AP84761



Package Size: #AP84761-1 50ug #AP84761-2 100ug #AP84761-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/threonine-protein kinase LMTK3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Apoptosis-associated tyrosine kinase 3, Lemur tyrosine kinase 3, Aatyk3
Accession No.	Uniprot ID:Q5XJV6
Uniprot	Q5XJV6
GeneID	381983;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Lmtk3

Note: This product is for in vitro research use only