

Lysyl oxidase homolog 2

Catalog No: #AP84762



Package Size: #AP84762-1 50ug #AP84762-2 100ug #AP84762-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysyl oxidase homolog 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Lysyl oxidase-like protein 2
Accession No.	Uniprot ID:B5DF27
Uniprot	B5DF27
GeneID	290350;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Loxl2

Note: This product is for in vitro research use only