

Epididymal-specific lipocalin-12

Catalog No: #AP84764



Package Size: #AP84764-1 50ug #AP84764-2 100ug #AP84764-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Epididymal-specific lipocalin-12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q6JVE5
Uniprot	Q6JVE5
GeneID	286256;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LCN12

Note: This product is for in vitro research use only