

Ligand-dependent nuclear receptor corepressor-like protein

Catalog No: #AP84765

Package Size: #AP84765-1 50ug #AP84765-2 100ug #AP84765-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ligand-dependent nuclear receptor corepressor-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LCoR-like protein,Mblk1-related protein 1,Transcription factor MLR1,Mr1
Accession No.	Uniprot ID:Q3U285
Uniprot	Q3U285
GeneID	209707;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lcor1

Note: This product is for in vitro research use only