

Interleukin-17A

Catalog No: #AP84767



Package Size: #AP84767-1 50ug #AP84767-2 100ug #AP84767-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--------------------------------------------------------|
| Product Name | Interleukin-17A |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | IL-17 |
| Accession No. | Uniprot ID:Q687Y7 |
| Uniprot | Q687Y7 |
| GeneID | 282863; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IL17A

Note: This product is for in vitro research use only