

Inhibin beta B chain

Catalog No: #AP84769



Package Size: #AP84769-1 50ug #AP84769-2 100ug #AP84769-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Inhibin beta B chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Activin beta-B chain
Accession No.	Uniprot ID:Q04999
Uniprot	Q04999
GeneID	16324;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Inhbb

Note: This product is for in vitro research use only