

Interleukin-2

Catalog No: #AP84771



Package Size: #AP84771-1 50ug #AP84771-2 100ug #AP84771-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	IL-2,T-cell growth factor,TCGF
Accession No.	Uniprot ID:Q29416
Uniprot	Q29416
GeneID	403989;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL2

Note: This product is for in vitro research use only