

# Phosphoribosyl pyrophosphate synthase-associated protein 1

Catalog No: #AP84772

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP84772-1 50ug #AP84772-2 100ug #AP84772-3 1mg

## Description

Product Name	Phosphoribosyl pyrophosphate synthase-associated protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PRPP synthase-associated protein 1,39 kDa phosphoribosypyrophosphate synthase-associated protein,PAP39
Accession No.	Uniprot ID:Q63468
Uniprot	Q63468
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Prpsap1

Note: This product is for in vitro research use only