

Interleukin-20 receptor subunit alpha

Catalog No: #AP84773



Package Size: #AP84773-1 50ug #AP84773-2 100ug #AP84773-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-20 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-20 receptor subunit alpha, Cytokine receptor class-II member 8, Cytokine receptor family 2 member 8, CRF2-8, IL-20R1, ZcytoR7, UNQ681/PRO1315
Accession No.	Uniprot ID:Q9UHF4
Uniprot	Q9UHF4
GeneID	53832;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IL20RA

Note: This product is for in vitro research use only