

Type I iodothyronine deiodinase

Catalog No: #AP84777



Package Size: #AP84777-1 50ug #AP84777-2 100ug #AP84777-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Type I iodothyronine deiodinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	5DI,DIOI,Type 1 DI,Type-I 5'-deiodinase,Itdi1,Txd1
Accession No.	Uniprot ID:P24389
Uniprot	P24389
GeneID	25430;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Dio1

Note: This product is for in vitro research use only