

Kallikrein-6

Catalog No: #AP84780



Package Size: #AP84780-1 50ug #AP84780-2 100ug #AP84780-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Kallikrein-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Neurosin,Protease M,SP59,Serine protease 18,Serine protease 9,Zyme,PRSS18,PRSS9
Accession No.	Uniprot ID:Q92876
Uniprot	Q92876
GeneID	5653;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KLK6

Note: This product is for in vitro research use only