

Calcium/calmodulin-dependent protein kinase type II subunit alpha

Catalog No: #AP84781

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84781-1 50ug #AP84781-2 100ug #AP84781-3 1mg

Description

Product Name	Calcium/calmodulin-dependent protein kinase type II subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CaM kinase II subunit alpha,CAMKA,KIAA0968
Accession No.	Uniprot ID:Q9UQM7
Uniprot	Q9UQM7
GeneID	815;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CAMK2A

Note: This product is for in vitro research use only