

Galectin-9

Catalog No: #AP84784



Package Size: #AP84784-1 50ug #AP84784-2 100ug #AP84784-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galectin-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Gal-9
Accession No.	Uniprot ID:Q3MHZ8
Uniprot	Q3MHZ8
GeneID	510813;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LGALS9

Note: This product is for in vitro research use only