

Ig gamma-1 chain C region secreted form

Catalog No: #AP84785



Package Size: #AP84785-1 50ug #AP84785-2 100ug #AP84785-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ig gamma-1 chain C region secreted form
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Igh-4
Accession No.	Uniprot ID:P01869
Uniprot	P01869
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ighg1

Note: This product is for in vitro research use only