

Interleukin-17A

Catalog No: #AP84786



Package Size: #AP84786-1 50ug #AP84786-2 100ug #AP84786-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-17A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-17,Cytotoxic T-lymphocyte-associated antigen 8,CTLA-8,Ctla8,Il17
Accession No.	Uniprot ID:Q62386
Uniprot	Q62386
GeneID	16171;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Il17a

Note: This product is for in vitro research use only